

产品名称: GRB2 (19X8) Rabbit Monoclonal Antibody  
产品货号: AMRe11750



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | GRB2 (19X8) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,FC,IF-P                  |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 纯化方式 (Purification) | Affinity purification                                                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                      |
|------------------------|------------------------------------------------------|
| 基因名 (Gene Name)        | GRB2                                                 |
| 别名 (Alternative Names) | ASH; Grb2; Grb3 3; HT027; NCKAP2; Protein ash; SEM5; |
| 基因 ID (Gene ID)        | 2885.0                                               |
| 蛋白 ID (SwissProt ID)   | P62993.                                              |

## 产品应用 (Application)

|                          |                                                                              |
|--------------------------|------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:20-1:50,FC 1:20-1:50,IF-P 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | 25kDa                                                                        |

## 研究背景 (Background)

**产品名称: GRB2 (19X8) Rabbit Monoclonal Antibody**  
**产品货号: AMRe11750**




Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.  
Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of GRB2 expression in HEK293 cell lysate.

## 注意事项 (Note)

For research use only .