

产品名称: GSTM1 (14O4) Rabbit Monoclonal Antibody
产品货号: AMRe11831



产品概述 (Summary)

产品名称 (Production Name)	GSTM1 (14O4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTM1
别名 (Alternative Names)	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; MGC26563; GSTM1a-1a; GSTM1b-1b;
基因 ID (Gene ID)	2944.0
蛋白 ID (SwissProt ID)	P09488.A synthetic peptide of human Glutathione S Transferase mu

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)

产品名称: GSTM1 (14O4) Rabbit Monoclonal Antibody
产品货号: AMRe11831




Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Involved in the formation of glutathione conjugates of both prostaglandin A2 (PGA2) and prostaglandin J2 (PGJ2) (PubMed:9084911). Participates in the formation of novel hepxilin regioisomers (PubMed:21046276).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of Glutathione S Transferase mu 1 in Jurkat,CHO-K1,HeLa cell lysates using Glutathione S Transferase mu 1 antibody(1:500 diluted).

注意事项 (Note)

For research use only .