

产品名称: GTF2I (1A17) Rabbit Monoclonal Antibody
产品货号: AMRe11843



产品概述 (Summary)

产品名称 (Production Name)	GTF2I (1A17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2I
别名 (Alternative Names)	BAP135; BTKAP1; DIWS; Gtf2i; IB291; SPIN; WBS; WBSCR6;
基因 ID (Gene ID)	2969.0
蛋白 ID (SwissProt ID)	P78347.A synthetic peptide of human TFII I

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	112kDa

研究背景 (Background)

Interacts with the basal transcription machinery by coordinating the formation of a multiprotein complex at the C-FOS

产品名称: GTF2I (1A17) Rabbit Monoclonal Antibody
产品货号: AMRe11843




promoter, and linking specific signal responsive activator complexes. Interacts with the basal transcription machinery by coordinating the formation of a multiprotein complex at the C-FOS promoter, and linking specific signal responsive activator complexes. Promotes the formation of stable high-order complexes of SRF and PHOX1 and interacts cooperatively with PHOX1 to promote serum-inducible transcription of a reporter gene derived by the C-FOS serum response element (SRE). Acts as a coregulator for USF1 by binding independently two promoter elements, a pyrimidine-rich initiator (Inr) and an upstream E-box. Required for the formation of functional ARID3A DNA- binding complexes and for activation of immunoglobulin heavy-chain transcription upon B-lymphocyte activation.

研究领域 (Research Area)

Basal transcription factors;

图片 (Image Data)



Western blot detection of TFII I in A549,HL-60 using TFII I antibody(1:1000 diluted)

注意事项 (Note)

For research use only .