

产品名称: GTPBP4 (14J8) Rabbit Monoclonal Antibody
产品货号: AMRe11848



产品概述 (Summary)

产品名称 (Production Name)	GTPBP4 (14J8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTPBP4
别名 (Alternative Names)	CRFG; GTP binding protein 4; Gtpbp4; NGB; NOG1;
基因 ID (Gene ID)	23560.0
蛋白 ID (SwissProt ID)	Q9BZE4.A synthetic peptide of human GTPBP4/NOG1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	74kDa

研究背景 (Background)

Involved in the biogenesis of the 60S ribosomal subunit. Involved in the biogenesis of the 60S ribosomal subunit.


产品名称: GTPBP4 (14J8) Rabbit Monoclonal Antibody
产品货号: AMRe11848



研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of GTPBP4/NOG1 in A549,HL-60 using GTPBP4/NOG1 antibody(1:1000 diluted)

注意事项 (Note)

For research use only .