

产品名称: Guanylyl Cyclase beta 1 (19P14) Rabbit Monoclonal Antibody
产品货号: AMRe11856



产品概述 (Summary)

产品名称 (Production Name)	Guanylyl Cyclase beta 1 (19P14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GUCY1B1
别名 (Alternative Names)	GCbeta1; GCS beta 1; GCS beta 3; GUC1B3; GUCSB3; GUCY1B1; Gucy1b1; GUCY1B3; SGC;
基因 ID (Gene ID)	2983.0
蛋白 ID (SwissProt ID)	Q02153.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	71kDa

**产品名称: Guanylyl Cyclase beta 1 (19P14) Rabbit
Monoclonal Antibody**
产品货号: AMRe11856




研究背景 (Background)

Mediates responses to nitric oxide (NO) by catalyzing the biosynthesis of the signaling molecule cGMP. Mediates responses to nitric oxide (NO) by catalyzing the biosynthesis of the signaling molecule cGMP.

研究领域 (Research Area)

Purine metabolism;Vascular smooth muscle contraction;Gap junction;Long-term depression;

图片 (Image Data)



Western blot analysis of Guanylyl Cyclase beta 1 expression in Human fetal brain lysate.

注意事项 (Note)

For research use only .