

产品名称: HAPLN1 (8N5) Rabbit Monoclonal Antibody
产品货号: AMRe11897



产品概述 (Summary)

产品名称 (Production Name)	HAPLN1 (8N5) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HAPLN1
别名 (Alternative Names)	Hyaluronan and proteoglycan link protein 1; Cartilage-linking protein 1; Cartilage-link protein; Proteoglycan link protein; HAPLN1; CRTL1
基因 ID (Gene ID)	1404.0
蛋白 ID (SwissProt ID)	P10915.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:200-1:2000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: HAPLN1 (8N5) Rabbit Monoclonal Antibody
产品货号: AMRe11897




研究背景 (Background)

HAPLN1 belongs to the HAPLN family. It contains 1 Ig like V type (immunoglobulin like) domain and 2 Link domains. This protein stabilizes the aggregates of proteoglycan monomers with hyaluronic acid in the extracellular cartilage matrix. Stabilizes the aggregates of proteoglycan monomers with hyaluronic acid in the extracellular cartilage matrix.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HAPLN1 expression in (1) Caco2 cell lysate; (2) Mouse spleen lysate.

注意事项 (Note)

For research use only .