

产品名称: Hemoglobin subunit alpha (15Z4) Rabbit Monoclonal Antibody
产品货号: AMRe11976



产品概述 (Summary)

产品名称 (Production Name)	Hemoglobin subunit alpha (15Z4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HBA1 HBA2; Alpha globin; Alpha-2 globin; HBH; Hemoglobin alpha chain;
别名 (Alternative Names)	Hemoglobin A1c; Hemoglobin subunit alpha; Hemoglobin alpha 2; Alpha-globin;
基因 ID (Gene ID)	3039;3040
蛋白 ID (SwissProt ID)	P69905.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,IF-P 1:200-1:500
----------------------	---

**产品名称: Hemoglobin subunit alpha (15Z4) Rabbit
Monoclonal Antibody**
产品货号: AMRe11976




蛋白分子量 (Molecular Weight) 15kDa

研究背景 (Background)

Involved in oxygen transport from the lung to the various peripheral tissues. Involved in oxygen transport from the lung to the various peripheral tissues.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Hemoglobin subunit alpha expression in K562 cell lysate.

注意事项 (Note)

For research use only .