

产品名称: HMGB1 (2P7) Rabbit Monoclonal Antibody
产品货号: AMRe12109



产品概述 (Summary)

产品名称 (Production Name)	HMGB1 (2P7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMGB1
别名 (Alternative Names)	
基因 ID (Gene ID)	3146.0
蛋白 ID (SwissProt ID)	P09429.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:20-1:100,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	25 kDa

研究背景 (Background)


产品名称: HMGB1 (2P7) Rabbit Monoclonal Antibody
产品货号: AMRe12109




研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of HMGB1 in K562, rat Brain, C6, 3T3, HeLa lysates using HMGB1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using HMGB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .