

产品名称: HNF1 beta (15Y3) Rabbit Monoclonal Antibody
产品货号: AMRe12125



产品概述 (Summary)

产品名称 (Production Name)	HNF1 beta (15Y3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNF1B
别名 (Alternative Names)	Hepatocyte nuclear factor 1 beta; HNF 1B; HNF 2; HNF-1-beta; HNF-1B; HNF1 beta; HNF1 homeobox B; HNF1B; HNF2; HPC11; LF B3; MODY 5; TCF 2 protein; VHNF 1;
基因 ID (Gene ID)	6928.0
蛋白 ID (SwissProt ID)	P35680.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:2000,IF-P 1:200-1:2000
蛋白分子量 (Molecular Weight)	61kDa


产品名称: HNF1 beta (15Y3) Rabbit Monoclonal Antibody
产品货号: AMRe12125

研究背景 (Background)

Belongs to the HNF1 homeobox family. Transcription factor, probably binds to the inverted palindrome 5'-GTAAATNATTAAC-3'. Transcription factor that binds to the inverted palindrome 5'-GTAAATNATTAAC-3' (PubMed:7900999, PubMed:17924661). Binds to the FPC element in the cAMP regulatory unit of the PLAU gene (By similarity). Transcriptional activity is increased by coactivator PCBD1 (PubMed:24204001).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HNF1 beta expression in Human kidney lysate.

注意事项 (Note)

For research use only .