

产品名称: HSD11B1 (17P16) Rabbit Monoclonal Antibody
产品货号: AMRe12219



产品概述 (Summary)

产品名称 (Production Name)	HSD11B1 (17P16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSD11B1
别名 (Alternative Names)	11-DH; 11-beta-HSD1; HSD11B1; HSD11; HSD11L; 11 beta HSD1;
基因 ID (Gene ID)	3290.0
蛋白 ID (SwissProt ID)	P28845.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:2000,ICC/IF 1:50-1:100,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	32kDa

研究背景 (Background)


产品名称: HSD11B1 (17P16) Rabbit Monoclonal Antibody
产品货号: AMRe12219



Catalyzes reversibly the conversion of cortisol to the inactive metabolite cortisone. Catalyzes reversibly the conversion of 7-ketocholesterol to 7-beta-hydroxycholesterol. In intact cells, the reaction runs only in one direction, from 7-ketocholesterol to 7-beta-hydroxycholesterol. Catalyzes reversibly the conversion of cortisol to the inactive metabolite cortisone (PubMed:10497248). Catalyzes reversibly the conversion of 7-ketocholesterol to 7-beta-hydroxycholesterol. In intact cells, the reaction runs only in one direction, from 7- ketocholesterol to 7-beta-hydroxycholesterol (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSD11B1 expression in human fetal liver lysate.

注意事项 (Note)

For research use only .