

产品名称: HSF2 (1X11) Rabbit Monoclonal Antibody
产品货号: AMRe12223



产品概述 (Summary)

产品名称 (Production Name)	HSF2 (1X11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSF2
别名 (Alternative Names)	Heat shock factor protein 2; HSF2; HSTF 2;
基因 ID (Gene ID)	3298.0
蛋白 ID (SwissProt ID)	Q03933.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: HSF2 (1X11) Rabbit Monoclonal Antibody
产品货号: AMRe12223



DNA-binding protein that specifically binds heat shock promoter elements (HSE) and activates transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked. DNA-binding protein that specifically binds heat shock promoter elements (HSE) and activates transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSF2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .