

产品名称: HSPA12A (5U3) Rabbit Monoclonal Antibody
产品货号: AMRe12264



产品概述 (Summary)

产品名称 (Production Name)	HSPA12A (5U3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA12A
别名 (Alternative Names)	FLJ13874 antibody;Heat shock 70 kDa protein 12A antibody;heat shock 70kD protein 12A antibody;heat shock 70kDa protein 12A antibody;HS12A_HUMAN antibody;Hspa12a antibody;KIAA0417 antibody
基因 ID (Gene ID)	259217.0
蛋白 ID (SwissProt ID)	O43301.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	75kDa

产品名称: HSPA12A (5U3) Rabbit Monoclonal Antibody
产品货号: AMRe12264




研究背景 (Background)

HSPA12A is an ATP binding protein. Adapter protein for SORL1, but not SORT1. Delays SORL1 internalization and affects SORL1 subcellular localization.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSPA12A expression in U-87 MG cell lysate.

注意事项 (Note)

For research use only .