

产品名称: Iba1 (9F7) Rabbit Monoclonal Antibody
产品货号: AMRe12328



产品概述 (Summary)

产品名称 (Production Name)	Iba1 (9F7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	IBA1; IRT1; AIF-1; IRT-1;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P55008/O70200.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:1000-1:10000,ICC/IF 1:100-1:200,FC 1:100-1:500
蛋白分子量 (Molecular Weight)	17kDa

研究背景 (Background)


产品名称: Iba1 (9F7) Rabbit Monoclonal Antibody
产品货号: AMRe12328



Actin-binding protein that enhances membrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play a role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Iba1 expression in Mouse testis cell lysate.

注意事项 (Note)

For research use only .