

产品名称: IDH1 (2P17) Rabbit Monoclonal Antibody
产品货号: AMRe12352



产品概述 (Summary)

产品名称 (Production Name)	IDH1 (2P17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IDH1
别名 (Alternative Names)	HEL-216; HEL-S-26; ICDH; IDCD; IDH; IDH1; IDP; IDPC; PICD;
基因 ID (Gene ID)	3417.0
蛋白 ID (SwissProt ID)	O75874.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:1000,FC 1:50-1:100,IP 1:20-1:50,IF-P 1:500-1:1000
蛋白分子量 (Molecular Weight)	47kDa

产品名称: IDH1 (2P17) Rabbit Monoclonal Antibody
产品货号: AMRe12352




研究背景 (Background)

Genetic variations are associated with cartilaginous tumors such as enchondroma or chondrosarcoma. Mutations of Arg-132 to Cys, Gly or His abolish the conversion of isocitrate to alpha-ketoglutarate. Instead, alpha-ketoglutarate is converted to R(-)-2-hydroxyglutarate.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of IDH1 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .