

产品名称: IgA (14V16) Rabbit Monoclonal Antibody
产品货号: AMRe12417



产品概述 (Summary)

产品名称 (Production Name)	IgA (14V16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	IGHA1; IGHA2; PIgR;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01876/P01877.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: IgA (14V16) Rabbit Monoclonal Antibody
产品货号: AMRe12417



Human IgA (immunoglobulin A) is a glycosylated protein of 160 kDa and is produced as a monomer or as a J chain linked dimer. Monomeric IgA constitutes 5-15 % of the serum immunoglobulins whereas dimeric IgA is localized to mucosa surfaces such as saliva, gastrointestinal secretion, bronchial fluids and milk. Mucosal IgA plays a major role in host defence by neutralising infectious agents at mucosal surfaces.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IgA expression in human plasma lysate.

注意事项 (Note)

For research use only .