

产品名称: IGFBP2 (12A1) Rabbit Monoclonal Antibody
产品货号: AMRe12429



产品概述 (Summary)

产品名称 (Production Name)	IGFBP2 (12A1) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFBP2
别名 (Alternative Names)	BP2; IBP2; IGF BP53; IGFBP2; IGFBP53;
基因 ID (Gene ID)	3485.0
蛋白 ID (SwissProt ID)	P18065.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:50-1:100,IP 1:20-1:50,IF-P 1:200-1:1000
蛋白分子量 (Molecular Weight)	35kDa

产品名称: IGFBP2 (12A1) Rabbit Monoclonal Antibody
产品货号: AMRe12429




研究背景 (Background)

Inhibits IGF-mediated growth and developmental rates. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors. Inhibits IGF-mediated growth and developmental rates. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of IGFBP2 expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .