

产品名称: IKK alpha/beta (10K2) Rabbit Monoclonal Antibody
产品货号: AMRe12472



产品概述 (Summary)

产品名称 (Production Name)	IKK alpha/beta (10K2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHUK/IKBKB
别名 (Alternative Names)	
基因 ID (Gene ID)	1147/3551
蛋白 ID (SwissProt ID)	O15111/O14920.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	85 kDa

研究背景 (Background)


产品名称: IKK alpha/beta (10K2) Rabbit Monoclonal Antibody
产品货号: AMRe12472



研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of IKK alpha/beta in 293T lysates using IKK alpha/beta antibody.

注意事项 (Note)

For research use only .