

产品名称: IL12 p40 (19Q5) Rabbit Monoclonal Antibody
产品货号: AMRe12490



产品概述 (Summary)

产品名称 (Production Name)	IL12 p40 (19Q5) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL12B
别名 (Alternative Names)	CLMF; CLMF p40; CLMF2; IL12B; NKSF; NKSF2;
基因 ID (Gene ID)	3593.0
蛋白 ID (SwissProt ID)	P29460.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)

产品名称: IL12 p40 (19Q5) Rabbit Monoclonal Antibody
产品货号: AMRe12490




Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated killer cells, and stimulate the production of IFN-gamma by resting PBMC. Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated killer cells, and stimulate the production of IFN-gamma by resting PBMC.

研究领域 (Research Area)

Cytokine-cytokine receptor interaction; Toll_Like; RIG-I-like receptor; Jak_STAT; Type I diabetes mellitus; Allograft rejection;

图片 (Image Data)



Western blot analysis of IL12 p40 expression in Human IL12 p40 recombinant protein.

注意事项 (Note)

For research use only .