

产品名称: IL13 (16E3) Rabbit Monoclonal Antibody
产品货号: AMRe12497



产品概述 (Summary)

产品名称 (Production Name)	IL13 (16E3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL13
别名 (Alternative Names)	ALRH; BHR1; IL13; Interleukin 13; NC30; P600;
基因 ID (Gene ID)	3596.0
蛋白 ID (SwissProt ID)	P35225.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	16kDa

研究背景 (Background)


产品名称: IL13 (16E3) Rabbit Monoclonal Antibody
产品货号: AMRe12497



Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327). Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327). May be critical in regulating inflammatory and immune responses (PubMed:8096327, PubMed:8097324). Positively regulates IL31RA expression in macrophages (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL13 expression in human recombinant protein cell lysate.

注意事项 (Note)

For research use only .