

产品名称: IL18 binding protein (18K17) Rabbit Monoclonal Antibody  
产品货号: AMRe12513



## 产品概述 (Summary)

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| 产品名称 (Production Name) | IL18 binding protein (18K17) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody                  |
| 宿主 (Host)              | Rabbit                                                  |
| 应用 (Application)       | WB                                                      |
| 种属反应性 (Reactivity)     | Human,Rat                                               |

## 产品性能 (Performance)

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 纯化方式 (Purification) | Affinity purification                                                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | IL18BP                                                                                                                               |
| 别名 (Alternative Names) | IL18 binding protein; IL18 Bpa; IL18BP; IL18BPα; Interleukin 18 binding protein; MC51L 53L 54L homolog gene product; Tadekinig alfa; |
| 基因 ID (Gene ID)        | 10068.0                                                                                                                              |
| 蛋白 ID (SwissProt ID)   | O95998.                                                                                                                              |

## 产品应用 (Application)

|                          |                  |
|--------------------------|------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000 |
| 蛋白分子量 (Molecular Weight) | 21kDa            |

**产品名称:** IL18 binding protein (18K17) Rabbit Monoclonal Antibody  
**产品货号:** AMRe12513




## 研究背景 (Background)

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response. Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of IL18 binding protein in HeLa cell lysate.

## 注意事项 (Note)

For research use only .