

产品名称: IL8 (6Z6) Rabbit Monoclonal Antibody
产品货号: AMRe12568



产品概述 (Summary)

产品名称 (Production Name)	IL8 (6Z6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CXCL8 CXCL8; GCP1; Granulocyte chemotactic protein 1; IL8; Interleukin 8; LECT; LUCT;
别名 (Alternative Names)	LYNAP; MDNCF; MONAP; NAF; NAP1; Neutrophil activating protein 1; SCYB8; TSG1;
基因 ID (Gene ID)	3576.0
蛋白 ID (SwissProt ID)	P10145.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	11kDa

产品名称: IL8 (6Z6) Rabbit Monoclonal Antibody
产品货号: AMRe12568




研究背景 (Background)

IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an inflammatory stimulus. IL-8(6-77) has a 5-10-fold higher activity on neutrophil activation, IL-8(5-77) has increased activity on neutrophil activation and IL-8(7-77) has a higher affinity to receptors CXCR1 and CXCR2 as compared to IL-8(1-77), respectively.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of IL8 expression in IL8 recombinant protein cell lysate.

注意事项 (Note)

For research use only .