

产品名称: Integrin alpha 9 (10J18) Rabbit Monoclonal Antibody
产品货号: AMRe12652



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha 9 (10J18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGA9
别名 (Alternative Names)	Alpha RLC; ITGA4L; ITGA9; RLC;
基因 ID (Gene ID)	3680.0
蛋白 ID (SwissProt ID)	Q13797 .

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	114kDa

产品名称: Integrin alpha 9 (10J18) Rabbit Monoclonal Antibody
产品货号: AMRe12652




研究背景 (Background)

Integrin alpha-9/beta-1 is a receptor for VCAM1, cytotactin and osteopontin. It recognizes the sequence A-E-I-D-G-I-E-L in cytotactin. Integrin alpha-9/beta-1 (ITGA9:ITGB1) is a receptor for VCAM1, cytotactin and osteopontin. It recognizes the sequence A-E-I-D- G-I-E-L in cytotactin.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Integrin alpha 9 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .