

产品名称: IRF2 (19Z15) Rabbit Monoclonal Antibody
产品货号: AMRe12738



产品概述 (Summary)

产品名称 (Production Name)	IRF2 (19Z15) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF2
别名 (Alternative Names)	Interferon regulatory factor 2; IRF 2;
基因 ID (Gene ID)	3660.0
蛋白 ID (SwissProt ID)	P14316.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	39kDa

研究背景 (Background)


产品名称: IRF2 (19Z15) Rabbit Monoclonal Antibody
产品货号: AMRe12738



Interferon regulatory factors (IRFs) comprise a family of transcription factors that function within the Jak/Stat pathway to regulate interferon (IFN) and IFN-inducible gene expression in response to viral infection. IRFs play an important role in pathogen defense, autoimmunity, lymphocyte development, cell growth, and susceptibility to transformation. IRF-2 acts as both a transcription activator and repressor. Specifically binds to the upstream regulatory region of type I IFN and IFN-inducible MHC class I genes (the interferon consensus sequence (ICS)) and represses those genes. Also acts as an activator for several genes including H4 and IL7. Constitutively binds to the ISRE promoter to activate IL7. Involved in cell cycle regulation through binding the site II (HiNF-M) promoter region of H4 and activating transcription during cell growth. Antagonizes IRF1 transcriptional activation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IRF2 expression in (1) HeLa cell lysate; (2) 3T3 cell lysate.

注意事项 (Note)

For research use only .