

产品名称: Lambda Light chain (14H15) Rabbit Monoclonal Antibody
产品货号: AMRe13188



产品概述 (Summary)

产品名称 (Production Name)	Lambda Light chain (14H15) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bence Jones protein; BJP; Constant region of lambda light chains; Ig lambda chain C regions; IGLC; IGLC1; IGLC2; IGLC3; IGLC6; IGLV; Immunoglobulin lambda constant 1; Mcg marker; Paraprotein;
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01701/P01842.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000,IHC 1:200-1:500,ICC/IF 1:100-1:200
----------------------	--

产品名称: Lambda Light chain (14H15) Rabbit Monoclonal Antibody
产品货号: AMRe13188




蛋白分子量 (Molecular Weight) 25kDa

研究背景 (Background)

All five immunoglobulin classes share the same basic four polypeptide chain structure of two heavy-chains and two light chains. There are five heavy chain types, and two light-chain types (Kappa and Lambda) both having a molecular weight of 22.5kDa. Any heavy-chain type can associate with either light-chain type, but on any immunoglobulin molecule both light-chains are of the same type.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Lambda Light chain expression in human plasma lysate.

注意事项 (Note)

For research use only .