

产品名称: Laminin gamma 1 (9A17) Rabbit Monoclonal Antibody
产品货号: AMRe13195



产品概述 (Summary)

产品名称 (Production Name)	Laminin gamma 1 (9A17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAMB2; lamc1; Laminin B2 chain; Laminin subunit gamma-1; Laminin-1 subunit gamma; Laminin-10 subunit gamma; Laminin-11 subunit gamma;
别名 (Alternative Names)	Laminin-2 subunit gamma; Laminin-3 subunit gamma; Laminin-4 subunit gamma; Laminin-6 subunit gamma; Laminin-7 subunit gamma; Laminin-8 subunit gamma; Laminin-9 subunit gamma;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P11047/Q5VYE7/Q6NVY8.Recombinant protein of human Laminin gamma 1

产品应用 (Application)

产品名称: Laminin gamma 1 (9A17) Rabbit Monoclonal Antibody
产品货号: AMRe13195




稀释比 (Dilution Ratio) WB 1:500-1:2000,IHC 1:200-1:500,FC 1:20-1:50,IP 1:20-1:50,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight) 178kDa

研究背景 (Background)

Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Laminin gamma 1 in HeLa, MCF-7, A549, HCT116 cell lysates using Laminin gamma 1 antibody (1:1000 diluted).

注意事项 (Note)

For research use only .