

产品名称: Liver Arginase (11G4) Rabbit Monoclonal Antibody
产品货号: AMRe13341



产品概述 (Summary)

产品名称 (Production Name)	Liver Arginase (11G4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Arg1
别名 (Alternative Names)	
基因 ID (Gene ID)	29221.0
蛋白 ID (SwissProt ID)	P07824.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000,IF-P 1:200-1:2000
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)


产品名称: Liver Arginase (11G4) Rabbit Monoclonal Antibody
产品货号: AMRe13341




研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Liver Arginase in rat liver, PC-12 lysates using Liver Arginase antibody.



Western blot analysis of Liver Arginase in mouse liver lysates using Liver Arginase antibody.

注意事项 (Note)

For research use only .