

产品名称: LMAN1 (16Y11) Rabbit Monoclonal Antibody
产品货号: AMRe13351



产品概述 (Summary)

产品名称 (Production Name)	LMAN1 (16Y11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMAN1
别名 (Alternative Names)	ERGIC53; F5F8D; FMFD1; Gp58; Lman1; MCFD1; MR60;
基因 ID (Gene ID)	3998.0
蛋白 ID (SwissProt ID)	P49257.A synthetic peptide of human LMAN1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:50-1:100,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	58kDa

研究背景 (Background)

Mannose-specific lectin. May recognize sugar residues of glycoproteins, glycolipids, or glycosylphosphatidyl inositol

产品名称: LMAN1 (16Y11) Rabbit Monoclonal Antibody
产品货号: AMRe13351




anchors and may be involved in the sorting or recycling of proteins, lipids, or both. The LMAN1-MCFD2 complex forms a specific cargo receptor for the ER-to-Golgi transport of selected proteins. Mannose-specific lectin. May recognize sugar residues of glycoproteins, glycolipids, or glycosylphosphatidyl inositol anchors and may be involved in the sorting or recycling of proteins, lipids, or both. The LMAN1-MCFD2 complex forms a specific cargo receptor for the ER-to-Golgi transport of selected proteins.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of LMAN1 in Hela,A549,HL-60,U251,U87-MG cell lysates using LMAN1 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .