

产品名称: LSP1 (8U19) Rabbit Monoclonal Antibody
产品货号: AMRe13465



产品概述 (Summary)

产品名称 (Production Name)	LSP1 (8U19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LSP1
别名 (Alternative Names)	Lymphocyte-specific protein 1; Protein pp52; 52 kDa phosphoprotein; Lymphocyte-specific antigen WP34; 47 kDa actin-binding protein;
基因 ID (Gene ID)	4046.0
蛋白 ID (SwissProt ID)	P33241.A synthetic peptide of human LSP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)


产品名称: LSP1 (8U19) Rabbit Monoclonal Antibody
产品货号: AMRe13465



LSP1 is a 52 kDa leukocyte-specific protein expressed in B cells, T cells, macrophages and neutrophils. The protein binds F-actin through its C-terminal actin-binding repeat and regulates cytoskeletal-related cellular movement such as polarization, chemotaxis and locomotion. LSP1 is a phospho-protein and has been identified as a major substrate of MAP kinase and PKC in lymphocytes and neutrophils. May play a role in mediating neutrophil activation and chemotaxis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Ramos cells using RM4626 at 1:1000.

注意事项 (Note)

For research use only .