

产品名称: MCL1 (18O13) Rabbit Monoclonal Antibody
产品货号: AMRe13712



产品概述 (Summary)

产品名称 (Production Name)	MCL1 (18O13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCL1
别名 (Alternative Names)	Bcl-2-like protein 3; Bcl-2-related protein EAT/mcl1; Bcl2-L-3; BCL2L3; EAT; Mcl-1; MCL1; mcl1/EAT; MCL1L; MCL1S;
基因 ID (Gene ID)	4170.0
蛋白 ID (SwissProt ID)	Q07820.A synthetic peptide of human MCL1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:100,ICC/IF 1:20-1:50,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	37kDa

研究背景 (Background)

产品名称: MCL1 (18O13) Rabbit Monoclonal Antibody
产品货号: AMRe13712




MCL1 a myeloid cell leukemia protein of the Bcl-2 family of proteins. Two alternatively spliced transcripts encoding distinct isoforms have been identified. The longer gene product (isoform 1) enhances cell survival by inhibiting apoptosis while the alternatively spliced shorter gene product (isoform 2) promotes apoptosis and is death-inducing. Involved in the regulation of apoptosis versus cell survival, and in the maintenance of viability but not of proliferation. Mediates its effects by interactions with a number of other regulators of apoptosis. Isoform 1 inhibits apoptosis. Isoform 2 promotes apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of MCL1 in C6 cell lysates using MCL1 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .