

产品名称: Mras (5X6) Rabbit Monoclonal Antibody
产品货号: AMRe14076



产品概述 (Summary)

产品名称 (Production Name)	Mras (5X6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRAS
别名 (Alternative Names)	Mras; R ras3; Ras related protein MRas; Ras related protein RRas3; RRas3; Xras;
基因 ID (Gene ID)	22808.0
蛋白 ID (SwissProt ID)	O14807.A synthetic peptide of human Mras

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)

May serve as an important signal transducer for a novel upstream stimuli in controlling cell proliferation. Weakly activates

产品名称: Mras (5X6) Rabbit Monoclonal Antibody
产品货号: AMRe14076




the MAP kinase pathway. Serves as an important signal transducer for a novel upstream stimuli in controlling cell proliferation. Activates the MAP kinase pathway.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of Mras in HeLa,Rat Brain lysates using Mras antibody(1:1000 diluted).

注意事项 (Note)

For research use only .