

产品名称: MYH11 (6C13) Rabbit Monoclonal Antibody
产品货号: AMRe14290



产品概述 (Summary)

产品名称 (Production Name)	MYH11 (6C13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH11 AAT4; FAA4; MYH11; Myosin 11; Myosin heavy chain 11; Myosin heavy chain smooth muscle isoform; Myosin heavy polypeptide 11 smooth muscle; SMHC; SMMHC;
别名 (Alternative Names)	
基因 ID (Gene ID)	4629.0
蛋白 ID (SwissProt ID)	P35749.A synthetic peptide of human Myosin , smooth muscle heavy chain 1 and 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:500,IF-P 1:100-1:500
蛋白分子量 (Molecular Weight)	227kDa

产品名称: MYH11 (6C13) Rabbit Monoclonal Antibody
产品货号: AMRe14290




研究背景 (Background)

Muscle contraction. Muscle contraction.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse spleen tissue using RM5298 at 1:1000.

注意事项 (Note)

For research use only .