

产品名称: MyoD1 (5M4) Rabbit Monoclonal Antibody
产品货号: AMRe14334



产品概述 (Summary)

产品名称 (Production Name)	MyoD1 (5M4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYOD1
别名 (Alternative Names)	Class C basic helix-loop-helix protein 1; MYF3; Myod 1; MYOD1; Myogenic differentiation 1; PUM;
基因 ID (Gene ID)	4654.0
蛋白 ID (SwissProt ID)	P15172.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	35kDa

产品名称: MyoD1 (5M4) Rabbit Monoclonal Antibody
产品货号: AMRe14334




研究背景 (Background)

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Activates muscle-specific promoters. Interacts with and is inhibited by the twist protein. This interaction probably involves the basic domains of both proteins. Acts as a transcriptional activator that promotes transcription of muscle-specific target genes and plays a role in muscle differentiation. Together with MYF5 and MYOG, co-occupies muscle-specific gene promoter core region during myogenesis. Induces fibroblasts to differentiate into myoblasts. Interacts with and is inhibited by the twist protein. This interaction probably involves the basic domains of both proteins (By similarity).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of MyoD1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .