

产品名称: Napsin A (14C2) Rabbit Monoclonal Antibody
产品货号: AMRe14406



产品概述 (Summary)

产品名称 (Production Name)	Napsin A (14C2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NAPSA
别名 (Alternative Names)	ASP4; KAP; Kdap; NAP1; NAPA; Napsa; Napsin-1; Napsin-A; Pronapsin A; SNAPA; TA01/TA02;
基因 ID (Gene ID)	9476.0
蛋白 ID (SwissProt ID)	O96009.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:500,IP 1:10-1:100,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight)	45kDa

产品名称: Napsin A (14C2) Rabbit Monoclonal Antibody
产品货号: AMRe14406




研究背景 (Background)

May be involved in processing of pneumocyte surfactant precursors. May be involved in processing of pneumocyte surfactant precursors.

研究领域 (Research Area)

Lysosome;

图片 (Image Data)



Western blot analysis of NAPSIN A expression in Human lung adenocarcinoma lysate.

注意事项 (Note)

For research use only .