

产品名称: NDRG1 (3L1) Rabbit Monoclonal Antibody
产品货号: AMRe14475



产品概述 (Summary)

产品名称 (Production Name)	NDRG1 (3L1) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDRG1
别名 (Alternative Names)	Protein NDRG1; Differentiation-related gene 1 protein; DRG-1; RTP; Rit42; NDRG1; CAP43; DRG1; RTP; targ1; TDD5; tdds;
基因 ID (Gene ID)	10397.0
蛋白 ID (SwissProt ID)	Q92597.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:50-1:200,ICC/IF 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	43kDa

产品名称: NDRG1 (3L1) Rabbit Monoclonal Antibody
产品货号: AMRe14475




研究背景 (Background)

May have a growth inhibitory role. Stress-responsive protein involved in hormone responses, cell growth, and differentiation. Acts as a tumor suppressor in many cell types. Necessary but not sufficient for p53/TP53-mediated caspase activation and apoptosis. Has a role in cell trafficking, notably of the Schwann cell, and is necessary for the maintenance and development of the peripheral nerve myelin sheath. Required for vesicular recycling of CDH1 and TF. May also function in lipid trafficking. Protects cells from spindle disruption damage. Functions in p53/TP53-dependent mitotic spindle checkpoint. Regulates microtubule dynamics and maintains euploidy.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of NDRG1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .