

产品名称: Nectin 2 (11O15) Rabbit Monoclonal Antibody
产品货号: AMRe14530



产品概述 (Summary)

产品名称 (Production Name)	Nectin 2 (11O15) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NECTIN2
别名 (Alternative Names)	HVEB; PRR2; CD112; PVRR2;
基因 ID (Gene ID)	5819.0
蛋白 ID (SwissProt ID)	Q92692.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	58kDa

研究背景 (Background)

产品名称: Nectin 2 (11O15) Rabbit Monoclonal Antibody
产品货号: AMRe14530




Probable cell adhesion protein. Modulator of T-cell signaling. Can be either a costimulator of T-cell function, or a coinhibitor, depending on the receptor it binds to. Upon binding to CD226, stimulates T-cell proliferation and cytokine production, including that of IL2, IL5, IL10, IL13, and IFNG. Upon interaction with PVRIG, inhibits T-cell proliferation. These interactions are competitive (PubMed:26755705). Probable cell adhesion protein (PubMed:9657005).

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Nectin 2 expression in SW480 cell lysate.

注意事项 (Note)

For research use only .