

产品名称: NEDD8 (5F5) Rabbit Monoclonal Antibody
产品货号: AMRe14538



产品概述 (Summary)

产品名称 (Production Name)	NEDD8 (5F5) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEDD8
别名 (Alternative Names)	NEDD8; Nedd-8;
基因 ID (Gene ID)	4738.0
蛋白 ID (SwissProt ID)	Q15843.A synthetic peptide of human NEDD8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:50-1:100,ICC/IF 1:20-1:50,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	9kDa

研究背景 (Background)

Neural precursor cell-expressed developmentally downregulated protein 8 (NEDD8), also known as Rub1 (related to

产品名称: NEDD8 (5F5) Rabbit Monoclonal Antibody
产品货号: AMRe14538




ubiquitin 1) in plants and yeast, is a member of the ubiquitin-like protein family. The covalent attachment of NEDD8 to target proteins, termed neddylation, is a reversible, multi-step process analogous to ubiquitination. Ubiquitin-like protein which plays an important role in cell cycle control and embryogenesis via its conjugation to a limited number of cellular proteins, such as cullins or p53/TP53 (PubMed:9694792, PubMed:10318914, PubMed:10597293, PubMed:11953428, PubMed:15242646, PubMed:14690597). Attachment of NEDD8 to cullins is critical for the recruitment of E2 to the cullin-RING-based E3 ubiquitin-protein ligase complex, thus facilitating polyubiquitination and proteasomal degradation of cyclins and other regulatory proteins (PubMed:9694792, PubMed:10318914, PubMed:10597293, PubMed:11953428, PubMed:20688984). Attachment of NEDD8 to p53/TP53 inhibits p53/TP53 transcriptional activity (PubMed:15242646). Covalent attachment to its substrates requires prior activation by the E1 complex UBE1C-APPBP1 and linkage to the E2 enzyme UBE2M (PubMed:14690597).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of NEDD8 in Jurkat cell lysates using NEDD8 antibody(1:1000 diluted).

产品名称: NEDD8 (5F5) Rabbit Monoclonal Antibody
产品货号: AMRe14538



注意事项 (Note)

For research use only .