

产品名称: NeuN (5W3) Rabbit Monoclonal Antibody
产品货号: AMRe14586



产品概述 (Summary)

产品名称 (Production Name)	NeuN (5W3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBFOX3
别名 (Alternative Names)	FLJ56884; FLJ58356; FOX3; HRNBP3;RBFOX3;NEUN; FOX-3;
基因 ID (Gene ID)	146713.0
蛋白 ID (SwissProt ID)	A6NFN3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:2000,IP 1:20-1:50,IF-P 1:500-1:2000
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

产品名称: NeuN (5W3) Rabbit Monoclonal Antibody
产品货号: AMRe14586




NeuN, also named as FOX3 and RBFOX3, is neuronal-specific nuclear protein with MW 38-50kd(~48kd). It is one of a family of 3 mammalian Fox homologues. FOX3(Or NeuN) functions in RNA-binding protein that regulates alternative splicing events. Pre-mRNA alternative splicing regulator. Regulates alternative splicing of RBFOX2 to enhance the production of mRNA species that are targeted for nonsense-mediated decay (NMD).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of NeuN expression in human fetal brain lysate.

注意事项 (Note)

For research use only .