

产品名称: NFE2 (14H9) Rabbit Monoclonal Antibody
产品货号: AMRe14645



产品概述 (Summary)

产品名称 (Production Name)	NFE2 (14H9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFE2
别名 (Alternative Names)	nfe2 ; Nuclear factor; p45 NFE2;
基因 ID (Gene ID)	4778.0
蛋白 ID (SwissProt ID)	Q16621.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:200-1:500,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)


产品名称: NFE2 (14H9) Rabbit Monoclonal Antibody
产品货号: AMRe14645



Component of the NF-E2 complex essential for regulating erythroid and megakaryocytic maturation and differentiation. Binds to the hypersensitive site 2 (HS2) of the beta-globin control region (LCR). Component of the NF-E2 complex essential for regulating erythroid and megakaryocytic maturation and differentiation. Binds to the hypersensitive site 2 (HS2) of the beta-globin control region (LCR). This subunit (NFE2) recognizes the TCAT/C sequence of the AP-1-like core palindrome present in a number of erythroid and megakaryocytic gene promoters. Requires MAFK or other small MAF proteins for binding to the NF-E2 motif. May play a role in all aspects of hemoglobin production from globin and heme synthesis to procurement of iron.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NFE2 expression in K562 cell lysate.

注意事项 (Note)

For research use only .