

产品名称: NFIB / NF1B2 (7L18) Rabbit Monoclonal Antibody
产品货号: AMRe14649



产品概述 (Summary)

产品名称 (Production Name)	NFIB / NF1B2 (7L18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFIB
别名 (Alternative Names)	CTF; HMGIC/NFIB; NF-I/B; NF1-B; NF1B; NF1B2; NFI-B; NFI-RED; Nfib; NFIB2; NFIB3;
基因 ID (Gene ID)	4781.0
蛋白 ID (SwissProt ID)	O00712.A synthetic peptide of human NFIB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	47kDa

产品名称: NFIB / NF1B2 (7L18) Rabbit Monoclonal Antibody
产品货号: AMRe14649


研究背景 (Background)

Recognizes and binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication. Transcriptional activator of GFAP, essential for proper brain development (PubMed:30388402). Recognizes and binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot detection of NFIB in C6, HeLa cell lysates using NFIB antibody (1:1000 diluted).

注意事项 (Note)

For research use only.