

产品名称: P2Y12 (3U6) Rabbit Monoclonal Antibody
产品货号: AMRe15608



产品概述 (Summary)

产品名称 (Production Name)	P2Y12 (3U6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2RY12 ADP glucose receptor; ADPGR; BDPLT8; HORK3; P2RY12; P2Y purinoceptor 12;
别名 (Alternative Names)	P2Y12; P2Y12 platelet ADP receptor; Platelet ADP receptor; Purinergic receptor P2RY12; Purinergic receptor P2Y12; SP1999;
基因 ID (Gene ID)	64805.0
蛋白 ID (SwissProt ID)	Q9H244.A synthetic peptide of human P2Y12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	40kDa

产品名称: P2Y12 (3U6) Rabbit Monoclonal Antibody
产品货号: AMRe15608




研究背景 (Background)

Receptor for ADP and ATP coupled to G-proteins that inhibit the adenylyl cyclase second messenger system. Not activated by UDP and UTP. Involved in platelets aggregation. Receptor for ADP and ATP coupled to G-proteins that inhibit the adenylyl cyclase second messenger system. Not activated by UDP and UTP. Required for normal platelet aggregation and blood coagulation.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot detection of P2Y12 in A549,HL-60,U251,U87-MG cell lysates using P2Y12 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .