

产品名称: P2Y6 (15H8) Rabbit Monoclonal Antibody
产品货号: AMRe15611



产品概述 (Summary)

产品名称 (Production Name)	P2Y6 (15H8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2RY6
别名 (Alternative Names)	P2ry6; P2Y6; P2Y6 receptor; PP2891;
基因 ID (Gene ID)	5031.0
蛋白 ID (SwissProt ID)	Q15077.A synthetic peptide of human P2Y6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:20-1:100,ICC/IF 1:20-1:50,IP 1:20-1:50,IF-P 1:20-1:100
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

产品名称: P2Y6 (15H8) Rabbit Monoclonal Antibody
产品货号: AMRe15611




Receptor for extracellular UDP > UTP > ATP. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system. Receptor for extracellular UDP > UTP > ATP. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of P2Y6 in K562, HeLa cell lysates using P2Y6 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.