

产品名称: pan Cadherin (16E7) Rabbit Monoclonal Antibody
产品货号: AMRe15721



产品概述 (Summary)

产品名称 (Production Name)	pan Cadherin (16E7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cadherin;CDH3;CDHP;Cadherin-1;Cadherin-2;Cadherin-3;Cadherin-4;Cadherin-5;CAM 120/80;CD144;CD324;CD325;Cdh4;CDH5;CDHE;CDHN;E-cadherin;N-cadherin;NCAD;Neural cadherin;P-cadherin;Uvomorulin;Vascular endothelial cadherin
别名 (Alternative Names)	

基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P12830/P19022/P22223/P33151/P55283/P39038.A synthetic peptide of human pan Cadherin

产品名称: pan Cadherin (16E7) Rabbit Monoclonal Antibody
产品货号: AMRe15721



产品应用 (Application)


稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	100kDa

研究背景 (Background)

Cadherins are members of a multigene family of single chain glycoprotein receptors mediating calcium dependent cell-cell adhesion. They play an important role in the growth and development of cells via the mechanisms of control of tissue architecture and the maintenance of tissue integrity. Cadherins are expressed in a tissue specific manner and are required for assembly of cells into solid tissue.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of pan Cadherin in HeLa,CHO-K1,C6,Rat Brain lysates using pan Cadherin antibody(1:500 diluted).

注意事项 (Note)

For research use only .