

产品名称: PAX8 (2K7) Rabbit Monoclonal Antibody
产品货号: AMRe15795



产品概述 (Summary)

产品名称 (Production Name)	PAX8 (2K7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAX8
别名 (Alternative Names)	Paired box 8; Paired box gene 8; paired box homeotic gene 8; Paired box protein Pax 8; Paired domain gene 8; PAX8;
基因 ID (Gene ID)	7849.0
蛋白 ID (SwissProt ID)	Q06710.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:500-1:1000
蛋白分子量 (Molecular Weight)	48kDa

产品名称: PAX8 (2K7) Rabbit Monoclonal Antibody
产品货号: AMRe15795




研究背景 (Background)

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells. Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

研究领域 (Research Area)

Pathways in cancer;Thyroid cancer;

图片 (Image Data)



Western blot analysis of PAX8 expression in SK-OV-3 cell lysate.

注意事项 (Note)

For research use only .