

产品名称: Plectin (16X18) Rabbit Monoclonal Antibody
产品货号: AMRe16263



产品概述 (Summary)

产品名称 (Production Name)	Plectin (16X18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLEC EBS1; EBSO; HD1; Hemidesmosomal protein 1; PCN; pleC; PLEC1 ;Plectin 1;
别名 (Alternative Names)	Plectin 6; Plectin; PLTN
基因 ID (Gene ID)	5339.0
蛋白 ID (SwissProt ID)	Q15149.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	532kDa

产品名称: Plectin (16X18) Rabbit Monoclonal Antibody
产品货号: AMRe16263




研究背景 (Background)

Interlinks intermediate filaments with microtubules and microfilaments and anchors intermediate filaments to desmosomes or hemidesmosomes. Could also bind muscle proteins such as actin to membrane complexes in muscle. May be involved not only in the filaments network, but also in the regulation of their dynamics. Interlinks intermediate filaments with microtubules and microfilaments and anchors intermediate filaments to desmosomes or hemidesmosomes. Could also bind muscle proteins such as actin to membrane complexes in muscle. May be involved not only in the filaments network, but also in the regulation of their dynamics. Structural component of muscle. Isoform 9 plays a major role in the maintenance of myofiber integrity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Plectin expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .