

产品名称: PMF1 (19K12) Rabbit Monoclonal Antibody
产品货号: AMRe16303



产品概述 (Summary)

产品名称 (Production Name)	PMF1 (19K12) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PMF1
别名 (Alternative Names)	Pmf1; Polyamine modulated factor 1;
基因 ID (Gene ID)	100527963;11243
蛋白 ID (SwissProt ID)	Q6P1K2.Recombinant protein of human PMF-1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	23kDa

研究背景 (Background)

Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore


产品名称: PMF1 (19K12) Rabbit Monoclonal Antibody
产品货号: AMRe16303



formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT. Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of PMF-1 in HeLa,A549 using PMF-1 antibody(1:1000 diluted)

注意事项 (Note)

For research use only .