

产品名称: Progesterone Receptor (4V13) Rabbit Monoclonal Antibody
产品货号: AMRe16517



产品概述 (Summary)

产品名称 (Production Name)	Progesterone Receptor (4V13) Rabbit Monoclonal Antibody
描述 (Description)	Rabbit Monoclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGR
别名 (Alternative Names)	NR3C3; PGR; PRGR; Progesterone receptor;
基因 ID (Gene ID)	5241.0
蛋白 ID (SwissProt ID)	P06401.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	99kDa

产品名称: Progesterone Receptor (4V13) Rabbit Monoclonal Antibody
产品货号: AMRe16517




研究背景 (Background)

Progesterone plays a central role in the reproductive events associated with the establishment and maintenance of pregnancy. Progesterone receptor, a member of the steroid receptor superfamily, mediates the physiologic effects of progesterone. The steroid hormones and their receptors are involved in the regulation of eukaryotic gene expression and affect cellular proliferation and differentiation in target tissues. Depending on the isoform, progesterone receptor functions as transcriptional activator or repressor.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Progesterone Receptor expression in T47 D cell lysate.

注意事项 (Note)

For research use only .