

产品名称: RAB22A (5R13) Rabbit Monoclonal Antibody
产品货号: AMRe16780



产品概述 (Summary)

产品名称 (Production Name)	RAB22A (5R13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB22A
别名 (Alternative Names)	Rab22; RAB22A;
基因 ID (Gene ID)	57403.0
蛋白 ID (SwissProt ID)	Q9UL26.A synthetic peptide of human RAB22A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:20-1:50, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	22kDa

研究背景 (Background)

Plays a role in endocytosis and intracellular protein transport. Mediates trafficking of TF from early endosomes to recycling

产品名称: RAB22A (5R13) Rabbit Monoclonal Antibody
产品货号: AMRe16780




endosomes (PubMed:16537905). Plays a role in endocytosis and intracellular protein transport. Mediates trafficking of TF from early endosomes to recycling endosomes (PubMed:16537905). Required for NGF-mediated endocytosis of NTRK1, and subsequent neurite outgrowth (PubMed:21849477). Binds GTP and GDP and has low GTPase activity. Alternates between a GTP-bound active form and a GDP-bound inactive form (PubMed:16537905).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of RAB22A in Rat Brain,C6,3T3,HeLa cell lysates using RAB22A antibody(1:1000 diluted).

注意事项 (Note)

For research use only .