

产品名称: RAC3 (9M14) Rabbit Monoclonal Antibody
产品货号: AMRe16824



产品概述 (Summary)

产品名称 (Production Name)	RAC3 (9M14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAC3
别名 (Alternative Names)	Rac1B; RAC3;
基因 ID (Gene ID)	5881.0
蛋白 ID (SwissProt ID)	P60763.A synthetic peptide of human RAC3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	21kDa

研究背景 (Background)

Plasma membrane-associated small GTPase which cycles between an active GTP-bound and inactive GDP-bound state.

产品名称: RAC3 (9M14) Rabbit Monoclonal Antibody
产品货号: AMRe16824




Plasma membrane-associated small GTPase which cycles between an active GTP-bound and inactive GDP-bound state. In active state binds to a variety of effector proteins to regulate cellular responses, such as cell spreading and the formation of actin-based protrusions including lamellipodia and membrane ruffles. Promotes cell adhesion and spreading on fibrinogen in a CIB1 and alpha-IIb/beta3 integrin-mediated manner.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of RAC3 in K562,Rat Brain cell lysates using RAC3 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .