

产品名称: RAD52 (13D10) Rabbit Monoclonal Antibody
产品货号: AMRe16847



产品概述 (Summary)

产品名称 (Production Name)	RAD52 (13D10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAD52
别名 (Alternative Names)	Rad52; RAD52 homolog;
基因 ID (Gene ID)	5893.0
蛋白 ID (SwissProt ID)	P43351.A synthetic peptide of human RAD52

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	46kDa

研究背景 (Background)

Involved in double-stranded break repair. Plays a central role in genetic recombination and DNA repair by promoting the

产品名称: RAD52 (13D10) Rabbit Monoclonal Antibody
产品货号: AMRe16847




annealing of complementary single-stranded DNA and by stimulation of the RAD51 recombinase. Involved in double-stranded break repair. Plays a central role in genetic recombination and DNA repair by promoting the annealing of complementary single-stranded DNA and by stimulation of the RAD51 recombinase.

研究领域 (Research Area)

Homologous recombination;

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using RM6417 at 1:1000.

注意事项 (Note)

For research use only .