

产品名称: RASA1 (10P8) Rabbit Monoclonal Antibody
产品货号: AMRe16907



产品概述 (Summary)

产品名称 (Production Name)	RASA1 (10P8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASA1
别名 (Alternative Names)	CMAVM; GAP; p120GAP; p120RASGAP; PKWS; RASA; RASA1; RasGAP;
基因 ID (Gene ID)	5921.0
蛋白 ID (SwissProt ID)	P20936.A synthetic peptide of human RASA1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:100,ICC/IF 1:20-1:50
蛋白分子量 (Molecular Weight)	116kDa

研究背景 (Background)

Inhibitory regulator of the Ras-cyclic AMP pathway. Stimulates the GTPase of normal but not oncogenic Ras p21. Inhibitory

产品名称: RASA1 (10P8) Rabbit Monoclonal Antibody
产品货号: AMRe16907




regulator of the Ras-cyclic AMP pathway. Stimulates the GTPase of normal but not oncogenic Ras p21; this stimulation may be further increased in the presence of NCK1.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of RASA1 in A549,HL-60,U87-MG cell lysates using RASA1 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .